RVR Dengue Combo NS1 & IgG/IgM Rapid Test (Serum/Plasma/Whole Blood)

REF RVR011
MDA Registration Number: IVDC7900739717

QUICK GUIDE


ASSAY PROCEDURE

- Step 1: Bring the specimen and test components to room temperature, if refrigerated or frozen. Mix the specimen well prior to assay once thawed.
- Step 2: When ready to test, open the pouch at the notch and remove device. Place the test device on a clean, flat surface.
- Step 3: Be sure to label the device with specimen's ID number.


FOR DETECTION OF DENGUE IgG/IgM

- Fill the capillary tube with specimen not to exceed the specimen line as shown in the images below.
- Holding the capillary tube vertically, dispense the entire specimen into the center of the sample well (**S well**) making sure that there are no air bubbles.
- Immediately add 3 drops (about 90-120 μL) of Sample Diluent into the buffer well (**B well**) with the bottle positioned vertically.

FOR DETECTION OF DENGUE NS1 ANTIGEN

- Fill the plastic dropper with specimen.
- Holding the dropper vertically, dispense 1 drop (about 30-45 μL) of serum/plasma or 2 drops (about 80-100 μL) of whole blood to the sample well (**S well**) making sure that there are no air bubbles.
- Immediately add 2 drops (about 80-100 μL) of Sample Diluent to the sample well (**S well**) with the bottle positioned vertically.

INTERPRETATION OF ASSAY RESULT

NEGATIVE RESULT	POSITIVE RESULTS	INVALID RESULTS
<p>If only the C line is present, the absence of any burgundy color in the G, M or T lines indicates that neither anti-dengue virus antibodies or dengue virus antigens are detected. The result is negative or non-reactive.</p>	<p>Samples with positive results should be confirmed with alternative testing method(s) and clinical findings before a positive determination is made.</p>	<p>If no C line is developed, the assay is invalid regardless of any burgundy color in the G, M or T lines as indicated below. Repeat the assay with a new device.</p>

MANUFACTURER

Chembio Diagnostics Malaysia Sdn. Bhd.

No. 2-1, 2nd Floor, N-TATT Building Jalan TP5, Taman Perindustrian UEP, Subang Jaya, 47600 Selangor Darul Ehsan Malaysia.
Tel: (603) 8084 4166 Email: myfeedback@chembio.com

© Copyright 2024, Chembio Diagnostics Malaysia Sdn Bhd. The Chembio logo is registered trademark of Chembio Diagnostic Systems, Inc.
Product Code: RVR011; Effective Date: 24 October 2024; Doc No.: RVR011-EN-MS-(QG) Rev 2

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




PENGELUAR

Chembio Diagnostics Malaysia Sdn. Bhd.

No. 2-1, 2nd Floor, N-TATT Building Jalan TP5, Taman Perindustrian UEP, Subang Jaya, 47600 Selangor Darul Ehsan Malaysia.

Tel: (603) 8084 4166 Emel: myfeedback@chembio.com

© Hak Cipta 2024, Chembio Diagnostics Malaysia Sdn Bhd. Logo Chembio adalah tanda dagangan berdaftar Chembio Diagnostic Systems, Inc.

Kod Produk: RVR011; Tarikh Berkuatkuasa: 24 Oktober 2024; No. Dokumen: RVR011-EN-MS-(QG) Rev 2