


FAST SCREENING of COVID-19
within 15-20 minutes!



RESZON COVID-19 RAPID ANTIBODY TEST

CE IVD FOR PROFESSIONAL USE ONLY, NOT FOR HOME TESTING



RESZON COVID-19 RAPID ANTIBODY TEST


FAST SCREENING of COVID-19 within 15-20 minutes

Test Description


This is an immunochromatographic assay designed for the qualitative detection of specific IgM and IgG antibodies against SARS-CoV-2, in human whole blood, serum or plasma. It is intended to be used by the professionals as a screening test and as an aid in the early diagnosis of SARS-CoV-2 infection.

Platform	Immunochromatographic Assay
Format	Cassette
Detection	Specific IgM and IgG antibodies against SARS-CoV-2
Specimen	Whole blood, serum or plasma
Specimen Volume	20 µl
Assay Time	15-20 minutes
Shelf Life	18 months

Procedures (Simplified Version)



Result Interpretation



Ordering Information

Cat. No	Product Description	Test Specimen	Product Format	Intended Use	Storage Temp	Packing Size
RCV-RD0101	Reszon COVID-19 Rapid IgG/IgM Test	Whole blood, serum or plasma	Cassette	Coronavirus disease COVID-19	4-30°C	25 T/ Kit
RCV-RD0101B	Reszon COVID-19 Rapid IgG/IgM Test (with disposable dropper, alcohol swab, lancet)	Whole blood, serum or plasma	Cassette	Coronavirus disease COVID-19	4-30°C	25 T/ Kit

COVID-19



Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)

A member of Revongen Corporation

Headquarters

Revongen Corporation Center, Level 17, Top Glove Tower, No. 16, Persiaran Setia Dagang, Setia Alam, Seksyen U13, 40170 Shah Alam, Selangor, Malaysia

+6 03 3359 1166

+6 03 3358 0303

info@reszonics.com

www.reszonics.com

Reszon Diagnostics

■ FOR EXPORT ONLY ■