

QUICK & EARLY detection of COVID-19 Infection



RESZON COVID-19 RAPID ANTIGEN TEST

CE IVD FOR PROFESSIONAL USE ONLY, NOT FOR HOME TESTING



RESZON COVID-19 RAPID ANTIGEN TEST


QUICK & EARLY detection of COVID-19 infection

Test Description


This is an immunochromatographic assay designed for the qualitative detection of specific SARS-CoV-2 antigen in human nasopharyngeal swab specimen. It is intended to be used by the professionals as a screening test and as an aid in the early diagnosis of SARS-CoV-2 infection.

| | |
|------------|-----------------------------|
| Platform | Immunochromatographic Assay |
| Format | Cassette |
| Detection | Specific SAR-CoV-2 antigen |
| Specimen | Nasopharyngeal swabs |
| Assay Time | 15-20 minutes |
| Shelf Life | 18 months |

Procedures (Simplified Version)



Result Interpretation



Ordering Information

| Cat. No | Product Description | Test Specimen | Product Format | Intended Use | Storage Temp | Packing Size |
|-------------|---|----------------------------|----------------|------------------------------|--------------|--------------|
| RCV-RD0201 | Reszon COVID-19 Rapid Antigen Test | Nasopharyngeal swab sample | Cassette | Coronavirus disease COVID-19 | 4-30°C | 25 T/ Kit |
| RCV-RD0201B | Reszon COVID-19 Rapid Antigen Test (with Nasopharyngeal Swab) | Nasopharyngeal swab sample | Cassette | Coronavirus disease COVID-19 | 4-30°C | 25 T/ Kit |

COVID-19



Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)

A member of Revongen Corporation

Headquarters

Revongen Corporation Center, Level 17, Top Glove Tower, No. 16, Persiaran Setia Dagang, Setia Alam, Seksyen U13, 40170 Shah Alam, Selangor, Malaysia

+6 03 3359 1166

+6 03 3358 0303

info@reszonics.com

www.reszonics.com

Reszon Diagnostics

■ FOR EXPORT ONLY ■