



THE FIRST MALAYSIAN MADE RAPID ANTIGEN HOME SELF-TESTING KIT

The non-intrusive COVID-19 rapid test kit for SELF-TEST
in human saliva and nasal swab samples.

Get your results at ease within 15 minutes from the
comfort of your home or any premise of your choice.

mySALIXIUM APP
CAN BE DOWNLOADED FROM



Available at
Nearest Pharmacies

**UJI SENDIRI,
LINDUNGI KOMUNITI**

#testyourselfprotectothers



Top Reasons To Choose SALIXIUM



AUTO-SYNC WITH MYSEJAHTERA

Every test result is automatically synchronized with MySejahtera app for further reporting to MOH.



TRACK & TRACE

The intelligently designed mySALIXIUM app allows every single test result to be tracked and traced accurately to ensure end-to-end traceability stemming from manufacturing site to end user and automatic reporting to Ministry of Health Malaysia.



FRIENDLY NOTIFICATION

mySALIXIUM app is smart and friendly to notify user the next course of action immediately upon completion of each test.



RESULT VALIDATION

Every single test cassette of SALIXIUM test kit is printed with a unique QR code for accurate and irreplicable result generation.

With the two-factor authentication technology using unique QR code and mySALIXIUM app, the end user unable to replicate or manipulate the test result using the same test cassette even without the presence of medical professional to validate the test result.



ENHANCE DETECTION RATE COMPARING TO JUST SALIVA SAMPLE ALONE

Self-collected saliva swab or nasal swab alone may miss certain detection of viral load. Combining both nasal and saliva swab samples have advantage to increase detection rate.



UNIQUE QR CODE AUTHENTICATION

Every single box of SALIXIUM test kit is affixed with distinctive "Halvec Verified" holographic QR code sticker that is printed with unique QR code to enable user to authenticate the product via mySALIXIUM app.

END TO END TRACEABILITY

1



2

DISTRIBUTOR



3

POINT OF PURCHASE



Support Local



Made in Malaysia



Endorsed by Ministry of Health, Malaysia



MDA Approved



Halvec Verified



Unique QR Code



mySALIXIUM App

SALIXIUM HIGHLIGHTS

SALIXIUM COVID-19 Rapid Antigen Test



No Pain

Less invasive sample collection method (saliva/nasal swab)



Fast

Rapid screening in less than 15 minutes



Unique QR Code

For accurate report generation, authentication and test tracing



Safe Swabs

Disposable swabs for saliva and nasal are included in individual sterile packing to minimize cross contamination



Waste Bag

Disposable waste bag is included for safe and easy disposal after use



Reliable

Endorsed by Ministry of Health, Malaysia



Halvec Verified

Tested by Halvec Laboratories, free from alcohol and animal DNA

4

CUSTOMER



5

END USER WHO USE THE SELF-TEST KIT



6

VISION AI FOR TEST VALIDATION



7

INTEGRATE WITH MYSEJAHTERA



8

REPORTING TO MOH





QUICK. ACCURATE. PAINLESS

**SCAN
ME**



**TUTORIAL VIDEO
QUICK USER GUIDIE**

Manufactured by

**Rezson Diagnostics
International Sdn. Bhd.**

 +6 03 3359 1166

 +6 03 3358 0303

 info@rezsonics.com

 www.rezsonics.com

PRODUCT SPECIFICATION

SALIXIUM Covid-19 Rapid Antigen Test



Sensitivity*

> 91%



Specificity*

100%



Test time

15 minutes



Storage condition

4-30°C



CE-IVD

Marked



Shelf life

18 months



Sample type

Saliva and anterior
nasal swab specimen

*vs FDA Emergency Use
Authorized real-time RT-PCR



Rezson Diagnostics

Contact Us

+6 016 216 0492

+6 016 828 8618