

RESZON

## COVID-19

RAPID ANTIGEN TEST

NASAL SWAB

REF RCV-RD0201N-01



SCAN ME



Please scan the QR code for more information



EASY

nasal self swab procedure



FAST

results in 15 minutes



RELIABLE

test performance

## PRODUCT SPECIFICATION

Sensitivity\*  
97.17%Test time  
15 minutesSample type  
Nasal swab specimenSpecificity\*  
100%Storage condition  
4-30°CShelf life  
18 months\* vs FDA Emergency Use  
Authorized real-time RT-PCR

For Self-testing

MADE IN MALAYSIA

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)  
A member of Revongen Corporation

## Headquarters

Revongen Corporation Center, Level 17, Top Glove Tower, No.16, Persiaran Setia Dagang, Setia Alam, Seksyen U13, 40170 Selangor, Malaysia.

+6 03 3359 1166

sales@reszonics.com

www.reszonics.com

■ FOR EXPORT ONLY ■



Reszon Diagnostics