

ENGLISH

RESZON



# COVID-19

## RAPID ANTIGEN TEST


SCAN ME



Please scan the QR code for more information

SALIVA

REF RCV-RD0201S-01



**EASY**  
saliva self  
collection



**FAST**  
results in  
15 minutes



**RELIABLE**  
test  
performance

### PRODUCT SPECIFICATION



**Sensitivity\***  
95.10%



**Test time**  
15 minutes



**Sample type**  
Saliva specimen



**Specificity\***  
100%



**Storage condition**  
4-30°C



**Shelf life**  
18 months

\* vs FDA Emergency Use  
Authorized real-time RT-PCR



For Self-testing

MADE IN MALAYSIA

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)  
A member of Revongen Corporation

Headquarters  
Revongen Corporation Center, Level 17, Top Glove Tower, No.16, Persiaran Setia Dagang,  
Setia Alam, Seksyen U13, 40170 Selangor, Malaysia.

+6 03 3359 1166

sales@reszonics.com

www.reszonics.com



Reszon Diagnostics