


RESZON

DOA RAPID

MDA Reg. No.: IVDB6168023-130652



Available 15
drugs in variety
of formats



No sample
preparation
required



Fast results
in 5 minutes



Easy result interpretation
as each drug type is clearly
labeled and color-coded



Urine Sample
Specimen


RESZON
DOA RAPID SINGLE DRUG TEST


100
TESTS
KIT

- Amphetamine (AMP)
- Buprenorphine (BUP)
- Benzodiazepines (BZO)
- Cocaine (COC)
- Erimin (ER5)
- Synthetic cannabinoids (K2/Spice)
- Ketamine (KET)
- Mitragynine/Kratom (MTG)
- Methamphetamine (MET)
- Morphine (MOP)
- Methadone (MTD)
- Ecstasy (MDMA)
- Oxycodone (OXY)
- Phencyclidine (PCP)
- Marijuana (THC)

TEST APPLICATIONS

- Pre-employment testing
- Scheduled testing
- Random testing
- Post-accident or for-cause testing
- Treatment follow up testing


RESZON
DOA RAPID DRUG COMBO TEST


25
TESTS
KIT

Able to detect from 2 to 10 different
drug combinations at the same time


RESZON
DOA RAPID DRUG CUP TEST


1 TEST
KIT


- Use a self-contained, one-step closed system with a secure lid to ensure total integrity of urine specimen
- Comes with integrated key activator (no dipping or pouring required)
- Able to detect from 4 to 12 different drug combinations at the same time

IVD


For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA


RESZON DOA RAPID SINGLE DRUG TEST




RESZON DOA RAPID DRUG COMBO TEST



RESZON DOA RAPID DRUG CUP TEST



RESULT INTERPRETATION



It is recommended that all positive results to be confirmed by using an appropriate analytical technique

IVD

For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA

RESZON DOA RAPID SINGLE DRUG TEST

Cat. No.	Product Description	Test Specimen	Product Format	Cut Off (ng/ml)	Storage Temp.	Packing Size
RDT-RD-AMP-100	Reszon DOA Rapid, Single Drug Test, AMP (Amphetamine)	Urine	Dipstick	1000	4-30°C	100 T/ Kit
RDT-RD-BUP-100	Reszon DOA Rapid, Single Drug Test, BUP (Buprenorphine)	Urine	Dipstick	10	4-30°C	100 T/ Kit
RDT-RD-BZO-100	Reszon DOA Rapid, Single Drug Test, BZO (Benzodiazepines)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-COC-100	Reszon DOA Rapid, Single Drug Test, COC (Cocaine)	Urine	Dipstick	150	4-30°C	100 T/ Kit
RDT-RD-ER5-100	Reszon DOA Rapid, Single Drug Test, ER5 (Erimin)	Urine	Dipstick	1000	4-30°C	100 T/ Kit
RDT-RD-K2-100	Reszon DOA Rapid, Single Drug Test, K2/Spice (Synthetic cannabinoids)	Urine	Dipstick	25	4-30°C	100 T/ Kit
RDT-RD-KET-100	Reszon DOA Rapid, Single Drug Test, KET (Ketamine)	Urine	Dipstick	1000	4-30°C	100 T/ Kit
RDT-RD-MDMA-100	Reszon DOA Rapid, Single Drug Test, MDMA (Ecstasy)	Urine	Dipstick	500	4-30°C	100 T/ Kit
RDT-RD-MET-100	Reszon DOA Rapid, Single Drug Test, MET (Methamphetamine)	Urine	Dipstick	1000	4-30°C	100 T/ Kit
RDT-RD-MOP-100	Reszon DOA Rapid, Single Drug Test, MOP (Morphine)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-MTD-100	Reszon DOA Rapid, Single Drug Test, MTD (Methadone)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-MTG-100	Reszon DOA Rapid, Single Drug Test, MTG (Mitragynine/Kratom)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-OXY-100	Reszon DOA Rapid, Single Drug Test, OXY (Oxycodone)	Urine	Dipstick	100	4-30°C	100 T/ Kit
RDT-RD-PCP-100	Reszon DOA Rapid, Single Drug Test, PCP (Phencyclidine)	Urine	Dipstick	25	4-30°C	100 T/ Kit
RDT-RD-THC-100	Reszon DOA Rapid, Single Drug Test, THC (Marijuana)	Urine	Dipstick	50	4-30°C	100 T/ Kit

RESZON DOA RAPID DRUG COMBO TEST

Cat. No.	Product Description	Test Specimen	Product Format	Storage Temp.	Packing Size
RDT-RD02,03,04,05-25	Reszon DOA Rapid, Combo 2-in-1 to 5-in-1	Urine	Dip card	4-30°C	25 T/ Kit
RDT-RD06,07,08,09,10-25	Reszon DOA Rapid, Combo 6-in-1 to 10-in-1	Urine	Twin card	4-30°C	25 T/ Kit

RESZON DOA RAPID DRUG CUP TEST

Cat. No.	Product Description	Test Specimen	Product Format	Storage Temp.	Packing Size
RDT-RDC04,05,06,07	Reszon DOA Rapid, Cup 4-in-1 to 7-in-1	Urine	Cup	4-30°C	1 Test
RDT-RDC10	Reszon DOA Rapid, Multiplex Cup 10-in-1	Urine	Cup	4-30°C	1 Test
RDT-RDC12	Reszon DOA Rapid, Multiplex Cup 12-in-1	Urine	Cup	4-30°C	1 Test



For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA