

# FILARIASIS IgG4 RAPID TEST KIT

REF RLF-RD0101-25

MDA Reg. No.: IVDC3603823-154220



Detect specific IgG4 antibodies to

- *Wuchereria bancrofti*
- *Brugia malayi*
- *Brugia timori*



## HASSLE-FREE

no special  
equipment  
needed



## EASY

sample preparation  
protocol and visual  
result interpretation



## RELIABLE

simple to  
perform

## PRODUCT SPECIFICATION



**Sensitivity**  
96.53%



**Test time**  
10-15 minutes



**Sample volume**  
20 µl



**Sample type**  
Human Whole Blood /  
Serum / Plasma



**Specificity**  
99.57%



**Storage  
Condition**  
4-30°C



**Kit**  
25 tests



**Shelf life**  
24 months

IVD

For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA