

# RESZON DENGUE

## RAPID TEST KIT

MDA Reg. No.: IVDC2463523-144368

SCAN ME



Scan for more information



### Reszon Dengue NS1 Ag Rapid Test

Detect dengue virus NS1 antigen in human serum, plasma, or whole blood.



### Reszon Dengue Rapid IgG/ IgM Test

Detect specific IgM and IgG antibodies to dengue virus in human serum, plasma, or whole blood.

Presumptive diagnosis of primary, secondary, recent/recurrent dengue infection.



### Reszon Dengue Rapid Combo Dengue NS1 & IgG/IgM Test

Detect Dengue NS1 Antigens and specific IgM and IgG antibodies to dengue virus in human serum, plasma, or whole blood.

Presumptive diagnosis of primary, secondary, recent, recurrent, past, early recovery phase dengue infection.



#### HASSLE-FREE

no special equipment needed



#### FAST

results in 20 minutes



#### EASY

sample preparation protocol and visual result interpretation



#### RELIABLE

test performance



#### MINIMAL

sample volume used

IVD

For professional *in vitro* diagnostic use only.

# RESZON

# DENGUE

## RAPID TEST KIT

MDA Reg. No.: IVDC2463523-144368



### PRODUCT SPECIFICATION

**RDG-RD0202-10**  
**RDG-RD0202-25**  
Reszon Dengue  
NS1 Ag Rapid Test

**RDG-RD0102-10**  
**RDG-RD0102-25**  
Reszon Dengue Rapid  
IgG/ IgM Test

**RDG-RD0301-10**  
**RDG-RD0301**  
Reszon Dengue Rapid  
Combo Dengue NS1 &  
IgG/IgM Test

**Sensitivity** 95.80%

96.72%

NS1: 95.80%  
IgG/IgM: 96.72%

**Specificity** 96.15%

97.58%

NS1: 96.15%  
IgG/IgM: 97.58%

**Storage condition** 4-30°C

4-30°C

4-30°C

**Specimen volume** 75µL

10µL

75 µL  
(Dengue NS1 Test)  
10 µL  
(Dengue IgG/IgM Test)

**Shelf life\*** 24 months

24 months

24 months

**Specimen** Human serum,  
plasma or  
whole blood

Human serum,  
plasma or  
whole blood

Human serum,  
plasma or  
whole blood

**Assay time** 15-20 mins

15-20 mins

15-20 mins

\* From date of manufacture

IVD

For professional *in vitro* diagnostic use only.