



Scan for more information

MDA Reg. No.: IVDB6168023-130652



Available 15 drugs in variety of formats



No sample preparation required



Fast results in 5 minutes



Easy result interpretation as each drug type is clearly labeled and color-coded



Urine Sample Specimen


RESZON
DOA RAPID SINGLE DRUG TEST

100 TESTS KIT

- Amphetamine (AMP)
- Buprenorphine (BUP)
- Benzodiazepines (BZO)
- Cocaine (COC)
- Erimin (ER5)
- Synthetic cannabinoids (K2/Spice)
- Ketamine (KET)
- Mitragynine/Kratom (MTG)
- Methamphetamine (MET)
- Morphine (MOP)
- Methadone (MTD)
- Ecstasy (MDMA)
- Oxycodone (OXY)
- Phencyclidine (PCP)
- Marijuana (THC)


RESZON
DOA RAPID DRUG COMBO TEST

25 TESTS KIT

Able to detect from 2 to 10 different drug combinations at the same time


RESZON
DOA RAPID DRUG CUP TEST

1 TEST KIT
TEST APPLICATIONS



- Pre-employment testing
- Scheduled testing
- Random testing
- Post-accident or for-cause testing
- Treatment follow up testing
- Use a self-contained, one-step closed system with a secure lid to ensure total integrity of urine specimen
- Comes with integrated key activator (no dipping or pouring required)
- Able to detect from 4 to 12 different drug combinations at the same time

IVD


 For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA



RESZON DOA RAPID SINGLE DRUG TEST



RESZON DOA RAPID DRUG COMBO TEST




RESZON DOA RAPID DRUG CUP TEST




RESULT INTERPRETATION



Positive



Negative



Invalid

It is recommended that all positive results to be confirmed by using an appropriate analytical technique

IVD

For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA

RESZON DOA RAPID SINGLE DRUG TEST

| Cat. No. | Product Description | Test Specimen | Product Format | Cut Off (ng/ml) | Storage Temp. | Packing Size |
|-----------------|---|---------------|----------------|-----------------|---------------|--------------|
| RDT-RD-AMP-100 | Reszon DOA Rapid, Single Drug Test, AMP (Amphetamine) | Urine | Dipstick | 1000 | 4-30°C | 100 T/ Kit |
| RDT-RD-BUP-100 | Reszon DOA Rapid, Single Drug Test, BUP (Buprenorphine) | Urine | Dipstick | 10 | 4-30°C | 100 T/ Kit |
| RDT-RD-BZO-100 | Reszon DOA Rapid, Single Drug Test, BZO (Benzodiazepines) | Urine | Dipstick | 300 | 4-30°C | 100 T/ Kit |
| RDT-RD-COC-100 | Reszon DOA Rapid, Single Drug Test, COC (Cocaine) | Urine | Dipstick | 150 | 4-30°C | 100 T/ Kit |
| RDT-RD-ER5-100 | Reszon DOA Rapid, Single Drug Test, ER5 (Erimin) | Urine | Dipstick | 300 | 4-30°C | 100 T/ Kit |
| RDT-RD-K2-100 | Reszon DOA Rapid, Single Drug Test, K2/Spice (Synthetic cannabinoids) | Urine | Dipstick | 25 | 4-30°C | 100 T/ Kit |
| RDT-RD-KET-100 | Reszon DOA Rapid, Single Drug Test, KET (Ketamine) | Urine | Dipstick | 1000 | 4-30°C | 100 T/ Kit |
| RDT-RD-MDMA-100 | Reszon DOA Rapid, Single Drug Test, MDMA (Ecstasy) | Urine | Dipstick | 500 | 4-30°C | 100 T/ Kit |
| RDT-RD-MET-100 | Reszon DOA Rapid, Single Drug Test, MET (Methamphetamine) | Urine | Dipstick | 1000 | 4-30°C | 100 T/ Kit |
| RDT-RD-MOP-100 | Reszon DOA Rapid, Single Drug Test, MOP (Morphine) | Urine | Dipstick | 300 | 4-30°C | 100 T/ Kit |
| RDT-RD-MTD-100 | Reszon DOA Rapid, Single Drug Test, MTD (Methadone) | Urine | Dipstick | 300 | 4-30°C | 100 T/ Kit |
| RDT-RD-MTG-100 | Reszon DOA Rapid, Single Drug Test, MTG (Mitragynine/Kratom) | Urine | Dipstick | 300 | 4-30°C | 100 T/ Kit |
| RDT-RD-OXY-100 | Reszon DOA Rapid, Single Drug Test, OXY (Oxycodone) | Urine | Dipstick | 100 | 4-30°C | 100 T/ Kit |
| RDT-RD-PCP-100 | Reszon DOA Rapid, Single Drug Test, PCP (Phencyclidine) | Urine | Dipstick | 25 | 4-30°C | 100 T/ Kit |
| RDT-RD-THC-100 | Reszon DOA Rapid, Single Drug Test, THC (Marijuana) | Urine | Dipstick | 50 | 4-30°C | 100 T/ Kit |

RESZON DOA RAPID DRUG COMBO TEST

| Cat. No. | Product Description | Test Specimen | Product Format | Storage Temp. | Packing Size |
|-------------------------|---|---------------|----------------|---------------|--------------|
| RDT-RD02,03,04,05-25 | Reszon DOA Rapid, Combo 2-in-1 to 5-in-1 | Urine | Dip card | 4-30°C | 25 T/ Kit |
| RDT-RD06,07,08,09,10-25 | Reszon DOA Rapid, Combo 6-in-1 to 10-in-1 | Urine | Twin card | 4-30°C | 25 T/ Kit |

RESZON DOA RAPID DRUG CUP TEST

| Cat. No. | Product Description | Test Specimen | Product Format | Storage Temp. | Packing Size |
|--------------------|---|---------------|----------------|---------------|--------------|
| RDT-RDC04,05,06,07 | Reszon DOA Rapid, Cup 4-in-1 to 7-in-1 | Urine | Cup | 4-30°C | 1 T/ Kit |
| RDT-RDC10 | Reszon DOA Rapid, Multiplex Cup 10-in-1 | Urine | Cup | 4-30°C | 1 T/ Kit |
| RDT-RDC12 | Reszon DOA Rapid, Multiplex Cup 12-in-1 | Urine | Cup | 4-30°C | 1 T/ Kit |



For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA