

FILARIASIS IgG4 RAPID TEST KIT

REF RLF-RD0101-25

MDA Reg. No.: IVDC3603823-154220

SCAN ME



Please scan
the QR code
for more
information



Detect specific IgG4 antibodies to

- *Wuchereria bancrofti*
- *Brugia malayi*
- *Brugia timori*



HASSLE-FREE

no special
equipment
needed



EASY

sample preparation
protocol and visual
result interpretation



RELIABLE

simple to
perform

PRODUCT SPECIFICATION



Sensitivity
96.53%



Test time
10-15 minutes



Sample volume
20 µl



Sample type
Human Whole Blood /
Serum / Plasma



Specificity
99.57%



Storage
4-30°C



Kit
25 tests



Shelf life
24 months

IVD

For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA