



# LYMPHATIC FILARIASIS RAPID TEST KIT BRUGIA RAPID

REF RLF-RD0201

MDA Reg. No.: IVDC8841920-38490

SCAN ME



Please scan  
the QR code  
for more  
information



Detect specific IgG4 antibodies to  
*Brugia spp.* (*Brugia malayi* and *Brugia timori*)



## HASSLE-FREE

no special  
equipment  
needed



## EASY

sample preparation  
protocol and visual  
result interpretation



## RELIABLE

test  
performance

## PRODUCT SPECIFICATION



**Sensitivity**  
94.75%



**Test time**  
15-25 minutes



**Sample volume**  
35  $\mu$ l



**Sample type**  
Human Whole Blood /  
Serum / Plasma



**Specificity**  
99.26%



**Storage**  
4-30°C



**Kit**  
25 tests



**Shelf life**  
18 months

IVD

For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA