



**TEST YOURSELF,
PROTECT OTHERS**

#testyourselfprotectothers

SALIXIUM COVID-19 RAPID ANTIGEN TEST (ORAL & NASAL SWAB SAMPLES)

REF SLXHB1-0621001

MDA Reg. No.: IVDC5125422-99355



THE FIRST MALAYSIAN MADE COVID-19 RAPID ANTIGEN HOME SELF-TESTING KIT

The non-intrusive COVID-19 rapid test kit for SELF-TEST in human oral and nasal swab samples.

Get your results at ease within 15 minutes from the comfort of your home or any premise of your choice.



SCAN ME

TUTORIAL VIDEO
QUICK USER GUIDE

PRODUCT SPECIFICATION



Sensitivity*
97.24%



Test time
15 minutes



Sample type
Oral and anterior
nasal swab specimen



Specificity*
100%



Storage condition
4-30°C



Shelf life
18 months

* vs FDA Emergency Use
Authorized real-time RT-PCR

SALIXIUM HIGHLIGHTS



Fast
Rapid screening in less
than 15 minutes



Dual Swabs
Disposable swabs for oral and
nasal are included in individual
sterile packing to minimize
cross contamination



No Pain
Less invasive sample collection
method (oral & nasal swab)



Biohazard Waste Bag
Disposable biohazard waste bag
is included for safe and easy
disposal after use




Halvec Verified
Tested by Halvec Laboratories,
free from alcohol and animal DNA

MANUFACTURER

Reszon Diagnostics International Sdn. Bhd.

+603 8022 1162 info@reszonics.com www.reszonics.com

Reszon Diagnostics



Doc. No.: MM-04-04-11-EN
Rev. No.: 01
Effective Date: 2024.01.22