



# LYMPHATIC FILARIASIS RAPID TEST KIT BRUGIA RAPID

REF RLF-RD0201

MDA Reg. No.: IVDC8841920-38490

SCAN ME

Please scan  
the QR code  
for more  
information

Detect specific IgG4 antibodies to  
*Brugia spp.* (*Brugia malayi* and *Brugia timori*)

**HASSLE-FREE**no special  
equipment  
needed**EASY**sample preparation  
protocol and visual  
result interpretation**RELIABLE**test  
performance**PRODUCT SPECIFICATION****Sensitivity**  
94.75%**Test time**  
15-25 minutes**Sample volume**  
35  $\mu$ l**Sample type**  
Human Whole Blood /  
Serum / Plasma**Specificity**  
99.26%**Storage**  
4-30°C**Kit**  
25 tests**Shelf life**  
24 months

IVD

For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA