



COVID-19

RAPID ANTIGEN TEST

ORAL & NASAL SWAB SAMPLES

REF RCV-RD0201SN-01

MDA Reg. No.: IVDC4716022-109725

SCAN ME



Please scan the QR code for more information



EASY
oral and nasal self
swab procedures



FAST
results in
15 minutes



RELIABLE
test
performance

PRODUCT SPECIFICATION



Sensitivity*
97.24%



Test time
15 minutes



Sample type
Oral and Nasal
swab specimen



Specificity*
100%



Storage condition
4-30°C



Shelf life
24 months

* vs FDA Emergency Use
Authorized real-time RT-PCR

IVD

For Self-testing

MADE IN MALAYSIA

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)

A member of Hextar Healthcare Berhad

+603 8022 1162

info@reszonics.com

www.reszonics.com

Reszon Diagnostics



Doc. No.: MM-01-04-11-EN

Rev. No.: 03

Effective Date: 2025.02.13