

COVID-19

RAPID ANTIGEN TEST

SALIVA

REF RCV-RD0201S-01

MDA Reg. No.: IVDC1451223-122236

SCAN ME



Please scan the QR code for more information



EASY
saliva self
collection



FAST
results in
15 minutes



RELIABLE
test
performance

PRODUCT SPECIFICATION



Sensitivity*
95.10%



Test time
15 minutes



Sample type
Saliva specimen



Specificity*
100%



Storage condition
4-30°C



Shelf life
24 months

* vs FDA Emergency Use
Authorized real-time RT-PCR

IVD

For Self-testing

MADE IN MALAYSIA