



RESZON

ENGLISH

DOA RAPID

MDA Reg. No.: IVDB10354125-196558



Available 24
drugs in variety
of formats



No sample
preparation
required



Fast results
in 5 minutes



Easy result interpretation
as each drug type is clearly
labeled and color-coded



Urine Sample
Specimen



Available with
adulterant tests



RESZON

DOA RAPID SINGLE DRUG TEST



100
TESTS
KIT

- 6-Monoacetylmorphine (6-MAM)
- Amphetamine (AMP)
- Barbiturates (BAR)
- Buprenorphine (BUP)
- Benzodiazepines (BZO)
- Carisoprodol (CAR/SOMA)
- Cocaine (COC)
- Erimin (ER5)
- Ethyl Glucuronide (ETG)
- Fentanyl (FYL)
- Synthetic cannabinoids (K2/Spice)
- Ketamine (KET)
- Mitragynine/Kratom (MTG)
- Methamphetamine (MET)
- Morphine (MOP)
- Methadone (MTD)
- Ecstasy (MDMA)
- Oxycodone (OXY)
- Phencyclidine (PCP)
- Propoxyphene (PPX)
- Psilocybin (PY)
- Marijuana (THC)
- Tricyclic Antidepressants (TCA)
- Tramadol (TML)

TEST APPLICATIONS

- Pre-employment testing
- Scheduled testing
- Random testing
- Post-accident or for-cause testing
- Treatment follow up testing



RESZON

DOA RAPID DRUG COMBO TEST



25
TESTS
KIT

Able to detect from 2 to 10 different
drug combinations at the same time



RESZON

DOA RAPID DRUG CUP TEST



1 TEST
KIT

- Use a self-contained, one-step closed system with a secure lid to ensure total integrity of urine specimen
- Comes with integrated key activator (no dipping or pouring required)
- Able to detect from 4 to 12 different drug combinations at the same time

IVD

For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)
A member of Hextar Healthcare Berhad

+603 8022 1162

info@reszonics.com

www.reszonics.com


Reszon Diagnostics




Doc. No.: MM-01-03-05-EN(1)
Rev. No.: 01
Effective Date: 2025.12.08

© 2025 Reszon Diagnostics International Sdn. Bhd. All rights reserved. RESZON® is the registered trademark of Reszon Diagnostics International Sdn. Bhd.


RESZON® DOA RAPID SINGLE DRUG TEST




RESZON® DOA RAPID DRUG COMBO TEST



RESZON® DOA RAPID DRUG CUP TEST



RESULT INTERPRETATION



It is recommended that all positive results to be confirmed by using an appropriate analytical technique

RESZON® DOA RAPID SINGLE DRUG TEST

Cat. No.	Product Description	Test Specimen	Product Format	Cut Off (ng/ml)	Storage Temp.	Packing Size
RDT-RD-AMP-100	RESZON® DOA Rapid, Single Drug Test, AMP (Amphetamine)	Urine	Dipstick	1000	4-30°C	100 T/ Kit
RDT-RD-BAR-100	RESZON® DOA Rapid, Single Drug Test, BAR (Barbiturates)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-BUP-100	RESZON® DOA Rapid, Single Drug Test, BUP (Buprenorphine)	Urine	Dipstick	10	4-30°C	100 T/ Kit
RDT-RD-BZO-100	RESZON® DOA Rapid, Single Drug Test, BZO (Benzodiazepines)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-CAR-100	RESZON® DOA Rapid, Single Drug Test, CAR (Carisoprodol / SOMA)	Urine	Dipstick	1000	4-30°C	100 T/ Kit
RDT-RD-COC-100	RESZON® DOA Rapid, Single Drug Test, COC (Cocaine)	Urine	Dipstick	150	4-30°C	100 T/ Kit
RDT-RD-ER5-100	RESZON® DOA Rapid, Single Drug Test, ER5 (Erimin)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-ETG-100	RESZON® DOA Rapid, Single Drug Test, ETG (Ethyl Glucuronide)	Urine	Dipstick	500	4-30°C	100 T/ Kit
RDT-RD-FYL-100	RESZON® DOA Rapid, Single Drug Test, FYL (Fentanyl)	Urine	Dipstick	20	4-30°C	100 T/ Kit
RDT-RD-K2-100	RESZON® DOA Rapid, Single Drug Test, K2/Spice (Synthetic cannabinoids)	Urine	Dipstick	25	4-30°C	100 T/ Kit
RDT-RD-KET-100	RESZON® DOA Rapid, Single Drug Test, KET (Ketamine)	Urine	Dipstick	1000	4-30°C	100 T/ Kit
RDT-RD-MAM-100	RESZON® DOA Rapid, Single Drug Test, 6-MAM (6-Monoacetylmorphine)	Urine	Dipstick	10	4-30°C	100 T/ Kit
RDT-RD-MDMA-100	RESZON® DOA Rapid, Single Drug Test, MDMA (Ecstasy)	Urine	Dipstick	500	4-30°C	100 T/ Kit
RDT-RD-MET-100	RESZON® DOA Rapid, Single Drug Test, MET (Methamphetamine)	Urine	Dipstick	1000	4-30°C	100 T/ Kit
RDT-RD-MOP-100	RESZON® DOA Rapid, Single Drug Test, MOP (Morphine)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-MTD-100	RESZON® DOA Rapid, Single Drug Test, MTD (Methadone)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-MTG-100	RESZON® DOA Rapid, Single Drug Test, MTG (Mitragynine/Kratom)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-OXY-100	RESZON® DOA Rapid, Single Drug Test, OXY (Oxycodone)	Urine	Dipstick	100	4-30°C	100 T/ Kit
RDT-RD-PCP-100	RESZON® DOA Rapid, Single Drug Test, PCP (Phencyclidine)	Urine	Dipstick	25	4-30°C	100 T/ Kit
RDT-RD-PPX-100	RESZON® DOA Rapid, Single Drug Test, PPX (Propoxyphene)	Urine	Dipstick	300	4-30°C	100 T/ Kit
RDT-RD-PY-100	RESZON® DOA Rapid, Single Drug Test, PY (Psilocybin)	Urine	Dipstick	500	4-30°C	100 T/ Kit
RDT-RD-TCA-100	RESZON® DOA Rapid, Single Drug Test, TCA (Tricyclic Antidepressants)	Urine	Dipstick	1000	4-30°C	100 T/ Kit
RDT-RD-THC-100	RESZON® DOA Rapid, Single Drug Test, THC (Marijuana)	Urine	Dipstick	50	4-30°C	100 T/ Kit
RDT-RD-TML-100	RESZON® DOA Rapid, Single Drug Test, TML (Tramadol)	Urine	Dipstick	100	4-30°C	100 T/ Kit

RESZON® DOA RAPID DRUG CUP TEST

Cat. No.	Product Description	Test Specimen	Product Format	Storage Temp.	Packing Size
RDT-RDC04,05,06,07	RESZON® DOA Rapid, Cup 4-in-1 to 7-in-1	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC10	RESZON® DOA Rapid, Multiplex Cup 10-in-1	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC12	RESZON® DOA Rapid, Multiplex Cup 12-in-1	Urine	Cup	4-30°C	1 T/ Kit



For professional *in vitro* diagnostic use only.


MADE IN MALAYSIA

RESZON® DOA RAPID DRUG COMBO TEST

Cat. No.	Product Description	Test Specimen	Product Format	Storage Temp.	Packing Size
RDT-RD02,03,04,05,06-25	RESZON® DOA Rapid, Combo 2-in-1 to 6-in-1	Urine	Dip card	4-30°C	25 T/ Kit
RDT-RD07,08,09,10-25	RESZON® DOA Rapid, Combo 7-in-1 to 10-in-1	Urine	Twin card	4-30°C	25 T/ Kit

RESZON® DOA RAPID ADULTERANT TEST

Cat. No.	Product Description	Test Specimen	Product Format	Storage Temp.	Packing Size
RDT-RD06-AD1-25	RESZON® DOA Rapid, Combo 6-in-1 with adulterant tests (CRE, OXI, pH)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD06-AD2-25	RESZON® DOA Rapid, Combo 6-in-1 with adulterant tests (SG, NIT, GLU)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD06-AD3-25	RESZON® DOA Rapid, Combo 6-in-1 with adulterant tests (CRE, OXI, pH, SG, NIT, GLU)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD07-AD1-25	RESZON® DOA Rapid, Combo 7-in-1 with adulterant tests (CRE, OXI, pH)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD07-AD2-25	RESZON® DOA Rapid, Combo 7-in-1 with adulterant tests (SG, NIT, GLU)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD07-AD3-25	RESZON® DOA Rapid, Combo 7-in-1 with adulterant tests (CRE, OXI, pH, SG, NIT, GLU)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD08-AD1-25	RESZON® DOA Rapid, Combo 8-in-1 with adulterant tests (CRE, OXI, pH)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD08-AD2-25	RESZON® DOA Rapid, Combo 8-in-1 with adulterant tests (SG, NIT, GLU)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD08-AD3-25	RESZON® DOA Rapid, Combo 8-in-1 with adulterant tests (CRE, OXI, pH, SG, NIT, GLU)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD09-AD1-25	RESZON® DOA Rapid, Combo 9-in-1 with adulterant tests (CRE, OXI, pH)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RD09-AD2-25	RESZON® DOA Rapid, Combo 9-in-1 with adulterant tests (SG, NIT, GLU)	Urine	Twin card	4-30°C	25 T / Kit
RDT-RDC04-AD1	RESZON® DOA Rapid, Cup 4-in-1, with adulterant tests (CRE, OXI, pH)	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC04-AD2	RESZON® DOA Rapid, Cup 4-in-1, with adulterant tests (SG, NIT, GLU)	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC04-AD3	RESZON® DOA Rapid, Cup 4-in-1, with adulterant tests (CRE, OXI, pH, SG, NIT, GLU)	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC05-AD1	RESZON® DOA Rapid, Cup 5-in-1, with adulterant tests (CRE, OXI, pH)	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC05-AD2	RESZON® DOA Rapid, Cup 5-in-1, with adulterant tests (SG, NIT, GLU)	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC05-AD3	RESZON® DOA Rapid, Cup 5-in-1, with adulterant tests (CRE, OXI, pH, SG, NIT, GLU)	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC06-AD1	RESZON® DOA Rapid, Cup 6-in-1, with adulterant tests (CRE, OXI, pH)	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC06-AD2	RESZON® DOA Rapid, Cup 6-in-1, with adulterant tests (SG, NIT, GLU)	Urine	Cup	4-30°C	1 T/ Kit
RDT-RDC06-AD3	RESZON® DOA Rapid, Cup 6-in-1, with adulterant tests (CRE, OXI, pH, SG, NIT, GLU)	Urine	Cup	4-30°C	1 T/ Kit



For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA