

ABC COMBO TEST

FLU A+B & COVID-19 AG

NASAL SWAB

REF RCV-RD0301-01

MDA Reg. No.: IVDC5790924-181335

SCAN ME



Please scan the QR code for more information



EASY
nasal swab
procedure



FAST
results in
20 minutes



RELIABLE
test
performance

PRODUCT SPECIFICATION



Sensitivity

Influenza A: 91.25%
Influenza B: 92.05%
COVID-19: 96.05%



Test time

20 minutes



Sample type

Nasal swab
specimen



Specificity

100%

* vs FDA Emergency Use
Authorized real-time
RT-PCR



Storage condition

4-30°C



Shelf life

24 months

IVD

For professional *in vitro* diagnostic use

MADE IN MALAYSIA