

RESZON COVID-19 RAPID ANTIGEN TEST

CALITAN HIDUNG

REF RCV-RD0201N-01

MDA Reg. No.: IVD05654022-103435

BAHASA MELAYU

IMBAS SAYA



Imbas kod QR
untuk
maklumat lanjut



MUDAH
prosedur swab
calitan hidung



CEPAT
keputusan dalam
masa 15 minit



DIPERCAYAI
prestasi
ujian

SPESIFIKASI PRODUK



Sensitiviti*
97.17%



Tempoh masa ujian
15 minit



Jenis sampel
Calitan hidung



Spesifisiti*
100%



Suhu penyimpanan
4-30°C



Jangka hayat
24 bulan

* Berbanding dengan FDA
Emergency Use Authorized
real-time RT-PCR

IVD

Untuk Ujian Kendiri

BUATAN MALAYSIA

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)
A member of Hextar Healthcare Berhad

+603 8022 1162

info@reszonics.com

www.reszonics.com

Reszon Diagnostics



No. Dokumen: MM-01-04-02-MS
No. Semakan: 04
Tarikh Kuat Kuasa: 2025.12.08

© 2025 Reszon Diagnostics International Sdn. Bhd. Hak cipta terpelihara. RESZON® ialah tanda dagangan berdaftar Reszon Diagnostics International Sdn. Bhd.