

COVID-19

RAPID ANTIGEN TEST

SALIVA

REF RCV-RD0201S-01

MDA Reg. No.: IVDC1451223-122236

SCAN ME



Please scan the QR code for more information



EASY
saliva self
collection



FAST
results in
15 minutes



RELIABLE
test
performance

PRODUCT SPECIFICATION



Sensitivity*
95.10%



Test time
15 minutes



Sample type
Saliva specimen



Specificity*
100%



Storage condition
4-30°C



Shelf life
24 months

* vs FDA Emergency Use
Authorized real-time RT-PCR

IVD

For Self-testing

MADE IN MALAYSIA

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)
A member of Hextar Healthcare Berhad

+603 8022 1162

info@reszonics.com

www.reszonics.com



Doc. No.: MM-01-04-12-EN
Rev. No.: 04
Effective Date: 2025.12.08

© 2025 Reszon Diagnostics International Sdn. Bhd. All rights reserved. RESZON® is the registered trademark of Reszon Diagnostics International Sdn. Bhd.