

RESZON DENGUE

RAPID TEST KIT

ENGLISH

SCAN ME



Scan for more
information

MDA Reg. No.: IVDC2463523-144368



RESZON® Dengue NS1 Ag Rapid Test

Detect dengue virus NS1 antigen in human serum, plasma, or whole blood.



RESZON® Dengue Rapid IgG/ IgM Test

Detect specific IgM and IgG antibodies to dengue virus in human serum, plasma, or whole blood.

Presumptive diagnosis of primary, secondary, recent/recurrent dengue infection.



RESZON® Dengue Rapid Combo Dengue NS1 & IgG/IgM Test

Detect Dengue NS1 Antigens and specific IgM and IgG antibodies to dengue virus in human serum, plasma, or whole blood.

Presumptive diagnosis of primary, secondary, recent, recurrent, past, early recovery phase dengue infection.



HASSLE-FREE

no special
equipment
needed



FAST

results in
20 minutes



EASY

sample preparation
protocol and visual
result interpretation



RELIABLE

test
performance



MINIMAL

sample
volume used

IVD

For professional *in vitro* diagnostic use only.

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)
A member of Hextar Healthcare Berhad

+603 8022 1162

info@reszonics.com

www.reszonics.com

Reszon Diagnostics



Doc. No.: MM-01-04-06-EN
Rev. No.: 03
Effective Date: 2025.12.08

© 2025 Reszon Diagnostics International Sdn. Bhd. All rights reserved. RESZON® is the registered trademark of Reszon Diagnostics International Sdn. Bhd.

RESZON DENGUE

ENGLISH

RAPID TEST KIT

MDA Reg. No.: IVDC2463523-144368



PRODUCT SPECIFICATION

RDG-RD0202-10
RDG-RD0202-25
RESZON® Dengue
NS1 Ag Rapid Test

RDG-RD0102-10
RDG-RD0102-25
RESZON® Dengue
Rapid IgG/ IgM Test

RDG-RD0301-10
RDG-RD0301
RESZON® Dengue
Rapid Combo Dengue
NS1 & IgG/IgM Test

Sensitivity 95.80%

96.72%

NS1: 95.80%
IgG/IgM: 96.72%

Specificity 96.15%

97.58%

NS1: 96.15%
IgG/IgM: 97.58%

**Storage
condition** 4-30°C

4-30°C

4-30°C

**Specimen
volume** 75µL

10µL

75 µL
(Dengue NS1 Test)
10 µL
(Dengue IgG/IgM Test)

Shelf life* 24 months

24 months

24 months

Specimen Human serum,
plasma or
whole blood

Human serum,
plasma or
whole blood

Human serum,
plasma or
whole blood

Assay time 15-20 mins

15-20 mins

15-20 mins

* From date of manufacture

IVD

For professional *in vitro* diagnostic use only.

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)
A member of Hextar Healthcare Berhad

+603 8022 1162

info@reszonics.com

www.reszonics.com

Reszon Diagnostics



Doc. No.: MM-01-04-06-EN
Rev. No.: 03
Effective Date: 2025.12.08

© 2025 Reszon Diagnostics International Sdn. Bhd. All rights reserved. RESZON® is the registered trademark of Reszon Diagnostics International Sdn. Bhd.