



HCV RAPID TEST

REF RHC-RD0101-25, RHC-RD0101A-25, RHC-RD0101B-25

MDA Reg. No.: IVDD10681325-211023

SCAN ME



Please scan the QR code for more information



Qualitative detection of antibodies to Hepatitis C Virus (HCV) in human whole blood, serum or plasma specimens



HASSLE-FREE

no special equipment needed



FAST

results in 20 minutes



EASY

sample preparation protocol and visual result interpretation



RELIABLE

test performance



MINIMAL

sample volume used

PRODUCT SPECIFICATION

- Sensitivity** 98.8%
- Test time** 15-20 minutes
- Sample volume** 50 µl (whole blood) 25 µl (serum /plasma)
- Sample type** Human Whole Blood / Serum / Plasma
- Specificity** 99.1%
- Storage** 4-30°C
- Kit** 25 tests
- Shelf life** 24 months



For professional *in vitro* diagnostic use only.

MADE IN MALAYSIA