

RESZON TYPHOID

RAPID TEST KIT

REF RTF-RD0101 / RTF-RD0202

MDA Reg. No.: IVDC85711391719

SCAN ME



Please scan the QR code for more information



RTF-RD0101
RESZON® Typhoid Rapid IgM

Detect specific IgM antibodies against a specific *Salmonella typhi* OMP antigen.



RTF-RD0202
RESZON® Typhoid Rapid IgG/IgM

Detect IgM and IgG antibodies against specific *Salmonella typhi* OMP antigen.



HASSLE-FREE

no any special equipment needed



FAST

results in 20 minutes



EASY

sample preparation protocol and visual result interpretation



RELIABLE

test performance



MINIMAL

sample volume used

IVD

For professional *in vitro* diagnostic use only.

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)
A member of Hextar Healthcare Berhad

+603 8022 1162

info@reszonics.com

www.reszonics.com

Reszon Diagnostics



Doc. No.: MM-01-04-03-EN

Rev. No.: 04

Effective Date: 2025.12.08

© 2025 Reszon Diagnostics International Sdn. Bhd. All rights reserved. RESZON® is the registered trademark of Reszon Diagnostics International Sdn. Bhd.



RESZON

TYPHOID

RAPID TEST KIT

REF RTF-RD0101 / RTF-RD0202

MDA Reg. No.: IVDC85711391719



PRODUCT SPECIFICATION

	RTF-RD0101 RESZON® Typhoid Rapid IgM	RTF-RD0202 RESZON® Typhoid Rapid IgG/IgM
Sensitivity	85.9%	IgG: 94.44% IgM: 91.67%
Specificity	96.7%	IgG: 93.40% IgM: 95.24%
Storage condition	4-30°C	4-30°C
Specimen volume	35uL	45uL
Shelf life	24 months	24 months
Specimen	Human serum, plasma or whole blood	Human serum, plasma or whole blood
Assay time	15-20 mins (whole blood) 10-15 mins (serum/ plasma)	15-20 mins (serum/ plasma/ whole blood)

IVD

For professional *in vitro* diagnostic use only.

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)
A member of Hexstar Healthcare Berhad

+603 8022 1162

info@reszonics.com

www.reszonics.com



Doc. No.: MM-01-04-03-EN
Rev. No.: 04
Effective Date: 2025.12.08