




RELIABLE RAPID TESTING FOR CLINICIANS

RESPIRATORY . FILARIASIS .
HEPATITIS . DENGUE . HIV

Because Every Minute MATTERS

Reson
Diagnostics
INTERNATIONAL



ABC COMBO TEST

FLU A+B & COVID-19 AG

NASAL SWAB

MDA Reg. No.: IVDC5790924-181335

SCAN ME



Please scan the QR code for more information

Clinically validated by
French Laboratory



RESPIRATORY INFECTION

RAPID TESTING SOLUTION



Three in one differential diagnosis

One step nasal swab process differentiates 3 respiratory infections



Rapid

Results in 15 minutes



Enhanced scalability

Suitable for high-throughput testing programs, enabling large scale deployment in clinical settings



Sensitivity*

- 91.25% (Influenza A)
- 92.05% (Influenza B)
- 96.05% (COVID-19)



Specificity*

- 100% (Influenza A, Influenza B and COVID-19)

*Refer to product insert for details

Cat. No.	Product Description	Test Specimen	Product Format	Intended Use	Storage Temp.	Packing Size
RCV-RD0301-01	RESZON® ABC Combo Test - Flu A+B & COVID-19 Ag (Nasal Swab)	Nasal Swab	Cassette	Influenza A, B and COVID-19	4-30°C	1 T/ Kit

LYMPHATIC FILARIASIS INFECTION

RAPID TESTING SOLUTION



FILARIASIS IgG4

RAPID TEST KIT

MDA Reg. No.: IVDC3603823-154220

SCAN ME



Please scan the QR code for more information



Comprehensive screening

Qualitative detection of *Wuchereria bancrofti*, *Brugia malayi* and/or *Brugia timori* infections



Rapid

Results in 15 minutes



Large scale screening

Ideal for screening foreign workers from diverse endemic countries



Trusted globally

Validated by international research teams*



Sensitivity*

96.53%



Specificity*

99.57%

*Refer to product insert for details

Cat. No.	Product Description	Test Specimen	Product Format	Intended Use	Storage Temp.	Packing Size
RLF-RD0101-25	RESZON® Filariasis IgG4 Rapid Test Kit	Whole Blood/ Serum/ Plasma	Cassette	IgG4 antibodies specific to <i>Wuchereria bancrofti</i> , <i>Brugia malayi</i> and/or <i>Brugia timori</i>	4-30°C	25 T/ Kit

HEPATITIS INFECTION
RAPID TESTING SOLUTION



Easy and accurate screening
Qualitative identification of HBsAg positive individuals



Rapid
Results within 15-20 minutes



Sensitivity*
99.3%



Specificity*
99.8%

*Refer to product insert for details

 **RESZON**

HBsAg RAPID TEST

MDA Reg. No.: IVDD4340025-214346

SCAN ME



Please scan the QR code for more information



Cat. No.	Product Description	Test Specimen	Product Format	Intended Use	Storage Temp.	Packing Size
RHB-RD0101-25	RESZON® HBsAg Rapid Test	Whole Blood/ Serum/ Plasma	Cassette	Detection of Hepatitis B surface antigen (HBsAg)	4-30°C	25 T/ Kit

RHB-RD0101A-25 includes 25pcs of disposable droppers.
RHB-RD0101B-25 includes 25pcs of disposable droppers, alcohol swabs and safety lancets.



Easy and accurate screening
Qualitative identification of HCV positive individuals



Rapid
Results within 15-20 minutes



Sensitivity*
98.8%



Specificity*
99.1%

*Refer to product insert for details

 **RESZON**

HCV RAPID TEST

MDA Reg. No.: IVDD10681325-211023

SCAN ME



Please scan the QR code for more information



Cat. No.	Product Description	Test Specimen	Product Format	Intended Use	Storage Temp.	Packing Size
RHC-RD0101-25	RESZON® HCV Rapid Test	Whole Blood/ Serum/ Plasma	Cassette	Detection of specific antibodies to HCV	4-30°C	25 T/ Kit

RHC-RD0101A-25 includes 25pcs of disposable droppers.
RHC-RD0101B-25 includes 25pcs of disposable droppers, alcohol swabs and safety lancets.

RVR Dengue Combo NS1-IgG/IgM Rapid Test

MDA Registration No.: IVDC7900739717



Full Spectrum Dengue Detection

- Early-stage primary infection: NS1 positive
- Late-stage primary infection: IgM and/or IgG positive
- Previous infection: IgG positive
- Acute secondary infection: NS1, IgG, and/or IgM positive



Rapid

Results in 20 minutes



Sensitivity*

96.00%



Specificity*

96.4%

*Refer to product insert for details

SCAN ME



Please scan the QR code for more information

Cat. No.	Product Description	Test Specimen	Product Format	Intended Use	Storage Temp.	Packing Size
RVR-011	Dengue Combo NS1-IgG/IgM Rapid Test	Whole Blood/Serum/Plasma	Cassette	Detection and differentiation of IgG anti-dengue virus, IgM anti-dengue virus and dengue antigen (Dengue NS1 antigen)	2-30°C	25 T/ Kit

HIV INFECTION

RAPID TESTING SOLUTION

HIV 1/2 STAT-PAK® Assay

MDA Reg. No.: IVDD2866999918



WHO Prequalified



FDA Approved



CE-Marked



Performance Matters

The only HIV test that exceeds WHO performance thresholds¹



Small sample volume

5µl



Sensitivity

100%



Specificity

100%²



Built-in procedural control

- 1 Kosack CS et al. Journal of the International AIDS Society 2017, 20:21345
- 2 Refer to product insert for details

SCAN ME



Please scan the QR code for more information

Cat. No.	Product Description	Test Specimen*	Product Format	Intended Use	Storage Temp.	Packing Size
60-9500 -O	HIV 1/2 STAT-PAK® Assay	Whole Blood/Serum/Plasma	Cassette	Detection of antibodies of HIV Type 1 and 2	8-30°C	20 T/ Kit

Distributed by

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)

A member of Hextar Healthcare Berhad

info@reszonics.com +603 8022 1162 +6012 262 8689

www.reszonics.com Reszon Diagnostics

Reszon
Diagnostics
INTERNATIONAL

Doc. No.: MM-06-03-24-EN

Rev. No.: 02

Effective Date: 2025.12.10