



LYMPHATIC FILARIASIS RAPID TEST KIT

BRUGIA RAPID

REF RLF-RD0201

MDA Reg. No.: IVDC8841920-38490

For professional *in vitro* diagnostic use only.

SCAN ME



Please scan
the QR code
for more
information



Detect specific IgG4 antibodies to
Brugia spp. (*Brugia malayi* and *Brugia timori*)



HASSLE-FREE

no special
equipment
needed



EASY

sample preparation
protocol and visual
result interpretation



RELIABLE

test
performance

PRODUCT SPECIFICATION



Sensitivity
94.75%



Test time
15-25 minutes



Sample volume
35 µl



Sample type
Human Whole Blood /
Serum / Plasma



Specificity
99.26%



Storage
2-30°C



Kit
25 tests




Shelf life
24 months


IVD

RESZON® LYMPHATIC FILARIASIS RAPID TEST KIT - BRUGIA RAPID


FOR WHOLE BLOOD SAMPLE



FOR SERUM AND PLASMA SAMPLE



RESULT INTERPRETATION



REAGENTS AND MATERIALS SUPPLIED

- RESZON® Lymphatic Filariasis Rapid Test - Brugia Rapid (individually sealed aluminium pouch with a desiccant)
- Instruction for Use (product insert)
- Assay buffer

ORDERING INFORMATION

RESZON® LYMPHATIC FILARIASIS RAPID TEST KIT - BRUGIA RAPID

Product Code	Description	Packing Size
RLF-RD0201	RESZON® Lymphatic Filariasis Rapid Test Kit - Brugia Rapid	25 tests / kit

Contact us today to enquire: +6012 815 8689 | info@reszonics.com

Reszon Diagnostics International Sdn. Bhd. 201001033978 (917901-W)

A member of Hextar Healthcare Berhad

©2026 Reszon Diagnostics International Sdn. Bhd. All rights reserved. RESZON® is the registered trademark of Reszon Diagnostics International Sdn. Bhd. Doc No: MM-01-04-09-EN; Rev. No.: 06; Effective Date: 2026.08.10
The content of this material is intended for Healthcare Professionals only, not for general public.

www.reszonics.com

Reszon
Diagnostics
INTERNATIONAL

Reszon Diagnostics

